



Produced in the EU for:



Important Notice

Must only be given to patients with proven PKU under strict medical supervision.

Suitable from 3 years of age.

Not for use as a sole source of nutrition.

Diet must be supplemented with natural protein, water and other nutrients in prescribed quantities to supply phenylalanine, fluid and general nutritional requirements. For enteral use only.

Storage

Store in a cool, dry place.

Pack Size

4 x 110 gram.

A Food for Special Medical Purposes



PKU Easy Microtabs

PKU Easy Microtabs are a phenylalanine free protein substitute tablet containing essential and non essential amino acids.

PKU Easy Microtabs are suitable to be used in the dietary management of PKU patients aged 3 years and above.

Dosage PKU Easy Microtabs

Dosage to be determined by the clinician or dietitian and is dependent on the age, body weight and medical condition of the patient.

PKU Easy Microtabs are small coated (slow release) tablets which are to be swallowed with liquid.

Can also be taken together with custard, jam or fruit puree.

Due to the special coating, they have no smell and no taste.

If the lid of the PKU Easy Microtabs bottle is filled with tablets, the content is equal to 10 grams of protein.

14 grams of tablets = approx. 10 grams of protein.

Nutritional information

Nutritional Information per 100 grams	
Energy	1678 kJ
	396 kcal
Protein equivalent	70.8 g
Amino acids	85.04 g
Fat	3.6 g
of which saturated fat	3.59 g
Carbohydrate	13 g
of which sugars	0 g
Fibre	3.7 g
Salt	1.0 g

Amino acids per 100 grams		
L-Alanine	3.07	g
L-Arginine	4.9	g
L-Aspartate	7.76	g
L-Carnitine	0.08	g
L-Cystine	2.01	g
Glycine	7.68	g
L-Glutamine	6.02	g
L-Histidine	3.07	g
L-Isoleucine	5.31	g
L-Leucine	8.27	g
L-Lysine	5.5	g
L-Methionine	1.42	g
L-Proline	5.54	g
L-Serine	3.42	g
L-Threonine	5.31	g
L-Tyrosine	7.78	g
L-Valine	6.13	g
Taurine	0.12	g
L-Tryptophan	1.65	g
Total	85.04	g