

GLORIOUS ENTERPRISES

PRODUCT SHEET

Products Name: Richnau-Holmgren Ear Speculam, 4mm

4mm
∅



Glorious Art nr. 49.0286.02
Net weight: 0.060 KG (approximate)
Material: Stainless steel AISI 304
Finish: SATIN
Sterility: Non-Sterile (Re-Usable)

Exposure Temperature: 134°C (274°F)
Exposure Time: 3 Minutes

CE Classification: Class-I
Storage: 16°C to 30°C (Store the instruments in a closed cabinet, protected from dust, moisture and temperature fluctuations)

Shelf life: 3 years
life-span: 20 times of use:- after each use, product need to sterilize for ready to use in next operation or next surgery.

Disposal: No special decomposition or disposal required for Surgical Instruments as they don't contain any toxic or hazardous material