

CORFLO

GASTROSTOMY

PEG Tubes and Jejunal
Extension Tubes

with ENFit™ Connectors



The ONLY Polyurethane Traction Removable PEG Tube

CORFLO® PEG FEATURES

- Unique internal bumper design featuring collapsing polyurethane foam
 - Collapses during insertion
 - Allows PEG placement for patients with strictures or narrowing of the oesophagus
- Traction Removable
 - Bumper collapses during removal
 - Endoscopic removal procedures not required avoiding invasive procedure and avoiding the need for contraindicated practice of 'cut & push' method.
 - Saves Endoscopy time, unnecessary admissions with associated costs
- Polyurethane Material
 - Inner lumen larger and stronger than silicone
- Locking Y-adaptor
 - Maximum security: no leaks or inadvertent disconnections
 - Dual Port maintains closed system during flushing and medication administration
- Unique tube fixation device holds tube securely at right angle without kinking
 - Prevents movement post op while stoma is healing.
- Flexible Offering
 - Pull or Push (Over-the-Wire) kit options
- Distinctive CORFLO YELLOW radiopaque stripe

CORFLO® ULTRA Jejunal Tube and Over-the-Wire Jejunal Tube

JEJUNAL EXTENSION TUBE FEATURES

- CORFLO with Grip-Tip and Over-the Wire Jejunal Extension Tubes are designed specifically for use with CORFLO PEG tubes
- Centimetre length markings for accurate adjustment and placement.
 - **Grip-Tip** aids endoscopic placement
 - Adjust-A-Sleeve allows sizing of tube length for each patient, while providing a secure non-leaking fit within PEG tube
 - Gastric access for aspiration, decompression and medication is maintained
 - Entire tube is radiopaque improving fluoroscopic detection during insertion
 - Polyurethane provides a large inner diameter, reducing risk of blockage
 - Anti-clog distal port ensures long term placement of tube

The new CORFLO clear proximal Y with ENFit™ connectors is designed to help prevent feeding tube misconnections and enhance safety.*

- Clear material allows for visual inspection of connection with administration sets and for cleanliness before, during and after use
- Soft material and ergonomic, smooth shape was designed with human factors engineering and customer input
- Purple caps with easy-grip ridges and threaded caps to provide secure connections. Purple helps identify the tubes for Enteral Use.



CORPAK
MEDSYSTEMS

Ask for the **YELLOW** Tube...

*ISO 80369-compliant. Multiple patents pending.

U.S. PATENTS D750,757, D750,758, D750,759, AND OTHER PATENTS PENDING

	PRODUCT CODE	DESCRIPTION	UNIT OF ISSUE
PEG Tubes	50-6012E1	12 Fr, polyurethane, traction removable, dual port	Box of 2
	50-6016E1	16 Fr, polyurethane, traction removable, dual port	Box of 2
	50-6020E1	20 Fr, polyurethane, traction removable, dual port	Box of 2
	50-4016E1	Over the Wire 16 Fr, polyurethane, traction removable, dual port	Box of 2
	50-4020E1	Over the Wire 20 Fr, polyurethane, traction removable, dual port	Box of 2
	50-6112	Adaptor Repair Kit 12 Fr	Box of 5
	50-6116	Adaptor Repair Kit 16 Fr	Box of 5
	50-6120	Adaptor Repair Kit 20 Fr	Box of 5

	PRODUCT CODE	DESCRIPTION	UNIT OF ISSUE
Jejunal Extension Tubes	CORFLO-ULTRA JEJUNAL TUBE WITH PEG TWOOMEY ADAPTOR FOR USE WITH THE 16 Fr CORFLO PEG TUBE		
	50-7366	6 Fr x 91 cm, polyurethane, radiopaque	Box of 5
	CORFLO-ULTRA JEJUNAL TUBES WITH PEG TWOOMEY ADAPTOR FOR USE WITH ALL 20 Fr CORFLO PEG TUBES		
	51-7368	8 Fr x 91 cm polyurethane, radiopaque	Box of 5
	50-7361	10 Fr x 91 cm, polyurethane, radiopaque	Box of 5
	CORFLO-ULTRA OVER-THE-WIRE JEJUNAL TUBES WITH PEG TWOOMEY ADAPTOR FOR USE WITH ALL 20 Fr CORFLO-MAX PEG TUBES		
	50-1431	10 Fr x 109 cm, Over the Wire	Box of 2



Innovation | Quality | Safety

TO ORDER

IN THE UK

0800 144 4480
Fax: 0800 151 3191
orders@corpakuk.com
CORPAK MedSystems UK
The Beehive, City Place
Gatwick RH6 0PA

NORTHERN IRELAND HOSPITAL ACCOUNTS:

028 9056 7135
Fax: 028 9070 5623
patricia.cartlidge@sangers.co.uk

NORTHERN IRELAND OTHER ACCOUNTS:

028 9040 3409
Fax: 028 9040 1958
bradburycustomerservices@sangers.co.uk

THE REPUBLIC OF IRELAND

01 468 8456
Fax: 01 466 0782
allpharservices@unipharm.ie
Allphar Services Limited
Orders Department
4045 Kingswood Road Citywest
Business Park
Co. Dublin

For more product information, please visit
www.corpak.com

CORPAK University Training
& Education Portal:
www.corpakU.com

