




GLORIOUS ENTERPRISES

PRODUCT SHEET

Products Name: Stille-Mod Plaster Shear 23cm



Glorious Art nr. 58.0102.23

Finish: Satin

Net weight: 210-215 gram (approximate)

Material: Stainless steel AISI 420

Exposure Temperature: 134°C (274°F)

Exposure Time: 10 Minutes



P.O. Box 63, Naverland 11, DK 2600 Glostrup, Denmark
Tel: +45 43 62 26 93, Fax: +45 43 62 26 11 , Email: glorious@glorious-enterprises.dk