

Dr. Fog® Anti-Fog

- Offered with our easy squeeze bottle and sponge set or as a ready-to-use sponge, which is dry to touch and alcohol free
- Radiopaque
- Sterile, intended for single-use, and not made with natural rubber latex

Part No.	Description	Qty/Box
DF-3100	Original Solution with Sponge	20
DF-3100B	Original Solution with Sponge	300
DF-3120	Alcohol-Free Dry Treated Sponge	40



Elvis® Anti-Fog

- Offered in both a squeeze bottle and pretreated sponge
- Sterile, intended for single-use, and not made with natural rubber latex



Part No.	Description	Qty/Box
400750	Anti-Fog, Pre-Treated Sponge	10
400751	Anti-Fog Bottle with Sponge	10

Surgical Headband

- Used to wick sweat away from the clinician
- Non-sterile and not made with natural rubber latex



Part No.	Description	Qty/Box
400450	Surgical Headband	10

Endoscopic Kittner

- Our blunt dissecting instrument is 38cm long with a 5mm stainless steel shaft
- Available in both single and dual tip
- Radiopaque tip
- The brushed finish allows for safe use during laser surgery
- Sterile, intended for single-use, and not made with natural rubber latex

Part No.	Description	Qty/Box
KT-9101	Single Tip	20
KT-9105	Dual Tip	20



Sterion® Cautery Tip Cleaner


- Assists in the cleaning and removing of material from cautery tips during surgical procedures
- Consists of a foam pad with abrasive material on one side and strong adhesive on the other for attachment to both disposable and non-disposable drapes
- Cleans both mono-polar and bi-polar probes
- Radiopaque
- Sterile, intended for single-use, and not made with natural rubber latex



Part No.	Description	Qty/Box
032001PBX	Cautery Tip Cleaner	35
TC-100	Cautery Tip Cleaner	100
TC-100A	Cautery Tip Cleaner	25
TC-100B	Cautery Tip Cleaner	240

Sterion® and Richard-Allan® Vessel Loops

- Used in occluding, retracting, and identifying arteries, veins, tendons and nerves in surgical procedures
- Extruded from medical-grade silicone, with the exception of the Richard-Allan® vessel loops, which are made of high performance, medical grade elastomer
- Radiopaque
- Packaged both sterile and non-sterile
- Intended for single-use and not made with natural rubber latex



Sterion® Vessel Loops

2 Loops / Pkg; 10 Pkgs / Box Double Pouched

Part No.	Color	mm	Inches	
			Width x Length	Width x Length
011001PBX	Red Mini	0.83 x 406	0.03 x 16	
011002PBX	Blue Mini	0.83 x 406	0.03 x 16	
011003PBX	Yellow Mini	0.83 x 406	0.03 x 16	
011004PBX	White Mini	0.83 x 406	0.03 x 16	
011011PBX	Red Maxi	1.0 x 406	0.04 x 16	
011012PBX	Blue Maxi	1.0 x 406	0.04 x 16	
011013PBX	Yellow Maxi	1.0 x 406	0.04 x 16	
011014PBX	White Maxi	1.0 x 406	0.04 x 16	
011021PBX	White Supermaxi	1.3 x 559	0.05 x 22	
011022PBX	Blue Supermaxi	1.3 x 559	0.05 x 22	
011032PBX	Blue Micro	0.41 x 406	0.019 x 16	
011301PBX	Red Mini	0.83 x 762	0.03 x 30	
011302PBX	Blue Mini	0.83 x 762	0.03 x 30	
011303PBX	Yellow Mini	0.83 x 762	0.03 x 30	
011304PBX	White Mini	0.83 x 762	0.03 x 30	
011311PBX	Red Maxi	1.0 x 762	0.04 x 30	
011312PBX	Blue Maxi	1.0 x 762	0.04 x 30	
011313PBX	Yellow Maxi	1.0 x 762	0.04 x 30	
011314PBX	White Maxi	1.0 x 762	0.04 x 30	

Richard-Allan® Vessel Loops



2 Loops / Pkg; 50 Pkgs / Box Banded and Pouched

Part No.	Color	mm	Inches	
			Width x Length	Width x Length
3901	Maxi Red	1.0 x 457	0.04 x 18	
3902	Maxi Blue	1.0 x 457	0.04 x 18	
3903	Maxi White	1.0 x 457	0.04 x 18	
3904	Maxi Yellow	1.0 x 457	0.04 x 18	
3905	Mini Red	.83 x 457	0.03 x 18	
3906	Mini Blue	.83 x 457	0.03 x 18	
3907	Mini White	.83 x 457	0.03 x 18	
3908	Mini Yellow	.83 x 457	0.03 x 18	

4 Loops / Pkg; 25 Pkgs / Box

Part No.	Color	mm	Inches	
			Width x Length	Width x Length
3910	Four Color Maxi Set	1.0 x 457	0.04 x 18	
3912	Four Color Mini Set	0.83 x 457	0.03 x 18	