

PRODUKTDATABLAD

Infusomat Space Line Enteral m/overgang(8250833SP)



Attributnamn	Attributvärde
Brand name	Infusomat®
Product name (HC)	Infusomat® Space Line
Product type	Administration set
Application	Enteral Nutrition
Type of administration set	Pressure
Container connector	ENPlus, Multiconnector 3in1
Längd	320.00 cm
Length from pump outlet	215 cm
Inner-ø tube	3.00 mm
Outer-ø tube	4.10 mm
Type of filter	No filter
Position of filter	-

Attributnamn	Attributvärde
drops/ml	20
Air ventilation	vented
Flow control device	Roller clamp
Injection access	ENFit Stopcock, needle-free
with AirStop	No
with PrimeStop	No
Number of back check valves	-
Number of SafSite®	-
Number of needle-free valves	-
Connection	ENLock, Transition adapter, ENFit female
Light-protected	No
Color of light protection	-
Pump compatibility (Type)	Infusomat® Space, Infusomat® fm/fmS
Length of pump segment	100.00 mm
max. tryck	2.0 bar
Material	PUR
Catalogue attribute	Enteral Nutrition
Type of packaging	Medipeel
Product (unpack.) contains PVC	No
Product (unpack.) cont. Latex	No
Product (unpack.) cont. DEHP	No
Packaging contains PVC	No
Packaging contains Latex	No
Packaging contains DEHP	No
Product contains PVC	No
Product contains Latex	No
Product contains DEHP	No
Product contains Nickel	No
Product contains BPA	No

Attributnamn	Attributvärde
Disposable / Reusable	disposable
Sterility	sterile
Type of Sterilization	E0
Type of Sterilizat. (allowed)	E0
UMDNS-Version	1.0
UMDNS	11677-Enteral Feeding Kits
Priming volume	28.400 ml
Length (Inch)	126
Type of Injection Access	ENFIT Stopcock
eCl@ss Code (Rel. 11)	34220605
eCl@ss Description (Rel. 11)	Inf. pump tub. (electr. power)