



### **Product description**

Order number	55.526
Product description	Tube, working volume: 5 ml, (LxØ): 75 x 12 mm, material: PP, round base, transparent, push cap, 500 piece(s)/StackPack

### **Product characteristics**

Graduation	no
Closure type	push cap
Base shape	round base
Volume of work	5 ml
Centrifugation max (RCF)	5000 x g

### **Size**

Diameter	12 mm
Length of product	75 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material Polypropylene (PP)  
Colour of product transparent

### **Purity & certification**

Satisfies the requirement sealed at 0.6 bar/10 min  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Batched yes

### **Packaging**

Minimum order qty. 4000  
Type of smallest StackPack  
subpackaging  
Piece(s) / inner box 500  
Piece (s) / outer case 4000  
Piece(s) / pallet 120000  
Depth of box 381 mm  
Width of box 77 mm  
Height of box 176 mm  
Depth of case 398 mm  
Width of case 358 mm  
Height of case 325 mm  
Case volume 0.0463 cbm  
Weight of product 0.0023 kg  
Weight of case 10.17 kg  
EAN of inner box 4038917082211  
EAN case 4038917001304

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

