

/ Tracheostomy Necktapes

Tube Fixations for tracheostomised patients

Order information

Tracheostomy Necktape with velcro fastener / For single use, material: velours coated foam, colour: blue

Patient	Length	REF	Box
Infant / Child	355 mm	40-0005-011	10
Child / Adult	475 mm	40-0005-033	10
Adult, extra large	525 mm	40-0005-044	10



Tracheostomy Necktape with velcro fastener / For single use, material: cotton coated foam, colour: skin coloured

Patient	Length	REF	Box
Infant / Child	340 mm	40-0008-011	10
Child / Adult	490 mm	40-0008-033	10
Adult, extra large	590 mm	40-0008-044	10



Tracheostomy Necktape with velcro fastener, separable / For single use, material: cotton coated foam, colour: skin coloured

Patient	Length	REF	Box
Child / Adult	400 mm	40-0066-033	10



/ Tracheostomy Dressing

- Breathable, acts as a second skin
- Absorbs moisture
- Avoids pressure sores

Order information

Tracheostomy Dressing / For single use, sterile

Size	REF	Box
9.0x9.0 cm	603207	10



Endotape

Oral / Nasal tracheal tube fixation

Order information

Endotape / With elastic, silicone coated ties and velcro fastener, for single use

Patient	Length	REF	Box
Child and Adult	540 mm	40-8018-044*	10

* Contains or presence of natural rubber latex



Tube Holder

Oral / Nasal tracheal tube fixation

Order information

Tube Holder / For single use

Patient	Length	REF	Box
Child and Adult	540 mm	40-8010-033	10
Adult	650 mm	40-8010-044	10



Biteblock

Oral tube fixation for tracheal tube, Laryngeal Tube and Laryngeal Mask

Tube clips into biteblock and is therefore protected from occlusion

Order information

Biteblock / For single use, sterile, material: EVA

Size	For Tracheal Tube	Size Laryngeal Tube	Size Laryngeal Mask	Biteblock	Biteblock with elastic necktape	Box
4	6.5 - 8.0 mm I.D.	-	2, 2.5	REF 40-60-040-1	REF 40-66-040-1	20
5	8.5 - 11.0 mm I.D.	2, 2.5	-	REF 40-60-050-1	REF 40-66-050-1	20
6	-	3, 4, 5	3, 4, 5	REF 40-60-060-1	REF 40-66-060-1	20

