

PROXIMATE® Skin Staplers

The number one choice in skin staplers*

The **PROXIMATE® RH** Skin Stapler features a 360° rotating head assembly for improved visibility and a ratchet mechanism for easier staple placement.

The **PROXIMATE® PX** Skin Stapler provides many of the same features as the PROXIMATE® RH, but in a fixed-head format.

The **PROXIMATE® PLUS MD** Skin Stapler is a high-value, low-cost skin stapler that permits multi-directional release in an ergonomic design.



PROXIMATE® Rotating-Head Skin Stapler

PROXIMATE Rotating-Head Skin Stapler

- Head rotates 360° for improved visibility and access
- Ratchet mechanism for easy staple placement
- Strong, rectangular stainless steel staples resist opening and minimize staple rotation
- Ergonomic pistol-grip handle is intuitive and comfortable to use

PROXIMATE Fixed-Head Skin Stapler

- Ratchet Mechanism for easy staple placement
- Ergonomic pistol-grip handle is intuitive and comfortable to use



PROXIMATE® Fixed-Head Skin Stapler

PROXIMATE PLUS Multi-Directional Release Skin Stapler





- Improved kick-off spring design for multi-directional release
- Ergonomic pistol-grip handle is intuitive and comfortable to use
- Alignment indicator improves visibility



PROXIMATE® PLUS MD Fixed-Head Skin Stapler

Ordering Information

Skin Staplers

Description	Device	Device Code	Number of Staples	Staple Leg Length	Staple Wire Diameter	Staple Crown	Qty/box
PROXIMATE® Rotating-Head Skin Stapler		PRR35	35	3.9mm	0.53mm	5.7mm	6
		PRW35	35	3.9mm	0.58mm	6.9mm	6
PROXIMATE® PX Fixed-Head Skin Stapler		PXR35	35	3.9mm	0.53mm	5.7mm	6
		PXW35	35	3.9mm	0.58mm	6.9mm	6
PROXIMATE® Plus MD Fixed-Head Stapler		PMR35	35	3.9mm	0.53mm	5.7mm	6
		PMW35	35	3.9mm	0.58mm	6.9mm	6
Description	Device	Device Code	Qty/box				
Squeeze-Handle Staple Extractor		PSX	12				

Squeeze-Handle Staple Extractor is compatible with all PROXIMATE Skin Staplers

How to Order

All purchase orders are made to Johnson & Johnson HealthCare Systems, Inc. (JJHCS) or your authorized distributor.

Electronic Ordering Options

(Note: Placing orders electronically avoids minimum order fees for your hospital.)

Johnson & Johnson Gateway:
www.jnjgateway.com/commerce

For questions about your order, visit the website or call 1-866-JNJ-GATE (565-4283).

Global Healthcare Exchange (GHX):
www.ghx.com

For questions about your order, visit the website or call 1-800-YOUR-GHX (968-7449).

Electronic Data Interchange (EDI):
JJHCS EDI Help Line: 1-800-262-2888

Nonelectronic/Manual Ordering Options

Call JJHCS at 1-800-255-2500 between 8:30 a.m. and 9:00 p.m. Eastern Standard Time or fax your order to 1-732-562-2212.

Customer support

For product use assistance, clinical guidelines, service and repair, emergency assistance, copy of a 510(k) clearance letter or complaints, please contact our Customer Support Center at customersupport@eesus.jnj.com or by calling **1-877-ETHICON** (384-4266). Our Support Center is staffed 24 hours a day, 7 days a week by qualified nurses to answer your product-related questions.

For complete product details, see Instructions for Use.