

HUMIDIFICATION THE NEED FOR HEAT AND MOISTURE

HME-ONLY CONFIGURATIONS


The adult-paediatric HME small is available with tethered CO₂ port cap for increased safety.

All products are latex free.

Individually packed, sterile, in boxes of 25.


This is a selection of the Medtronic DAR™ HME range.

HME for Tracheostomised Patients




353/5921
with 200 cm oxygen tube

HME for Tracheostomised Patients




353/19004

HME for Tracheostomised Patients




353S19046

Adult-Paediatric HME SMALL




353S19007
353S19007TC
with tethered cap

Adult-Paediatric HME SMALL



353/5426
with PVC catheter mount

Adult-Paediatric Foam HME with Integrated Catheter Mount



353P5908