

Lectroflex

Triple-layer extension tube (PVC/EVA/PE) for IV catheters with 1 male Luer-lock end and 1 female Luer-lock end. Diameter: 2,5 x 4,0 mm - Pressure tested at 10 kg/cm2 (145 PSI)

Available in 7 lengths : 13.5 cm (code 220.010), 25 cm (code 220.020), 50 cm (code 220.050), 70 cm (code 220.070), 100 cm (code 220.100), 150 cm (code 220.150) and 200 cm (code 220.200).



Specifications

Additional information

Extension tube					
Code	Int. Ø mm	Ext. Ø mm	Prim. vol. ml	Length cm	
220.010	2.5	4.0	1.0	13.5	
220.020	2.5	4.0	1.60	25	
220.050	2.5	4.0	2.80	50	
220.070	2.5	4.0	3.80	70	
220.100	2.5	4.0	5.2	100	
220.150	2.5	4.0	7.7	150	
220.200	2.5	4.0	10	200	

Specifications

Additional information

- ✓ Contains Latex: No
- ✓ Contains DEHP: No
- ✓ Contains products of animal or biological origin: No
- ✓ Pyrogen-free: Yes