



Laerdal
helping save lives



The BAG II Disposable Resuscitator

The BAG II Disposable Resuscitator is Laerdal's option for disposable ventilation. For lung ventilation on patients, without clean up and the risk of cross-contamination, The BAG II Disposable Resuscitator is the top choice. The BAG II has 3 different options for ventilation bag volume: adult, child or infant.

Single units

- 845111 The BAG II Resuscitator Adult w/mask 5
- 845141 The BAG II Resuscitator Adult w/mask 4
- 845151 The BAG II Resuscitator Adult mask 5 w/inflation port
- 845152 The BAG II Resuscitator Child mask 3 w/inflation port
- 845121 The BAG II Resuscitator Child w/mask 3
- 845131 The BAG II Resuscitator Infant w/mask 1

Multipacks (12-Pack)

- 845211 The BAG II Resuscitator Adult w/mask 5
- 845241 The BAG II Resuscitator Adult w/mask 4
- 845251 The BAG II Resuscitator Adult w/mask 5
w/inf. port - pkg of 12
- 845252 The BAG II Resuscitator Child mask 3
w/inf. port - pkg of 12
- 845221 The BAG II Resuscitator Child w/mask 3
- 845223 The BAG II Resuscitator Child w/mask 2
- 845231 The BAG II Resuscitator Infant w/mask 1

Accessories

- 845240 Laerdal Disposable PEEP Valve Assembly (qty. 10)
- 845255 Laerdal Disposable Mask 1 w/inflation port (qty.20)
- 845256 Laerdal Disposable Mask 2 w/inflation port (qty.20)
- 845250 Laerdal Disposable Mask 3 w/inflation port (qty.20)
- 845260 Laerdal Disposable Mask 4 w/inflation port (qty.20)
- 845270 Laerdal Disposable Mask 5 w/inflation port (qty.20)

Laerdal Medical AB

Värmdövägen 84
131 54 Nacka
Sweden

Telefon: +46 85 56 14 610
Epost: kontakt@laerdal.se