









## Superset® catheter mounts

The Superset range of catheter mounts have a 'concertina' design for total positioning control. This allows the mount to be set in any position required. Deadspace can be adjusted from 20ml to 40ml (code 3504: 15-36ml) and length from 70mm to 150mm (excluding connectors).

	Code	Box qty	Length	Sterile option code	Sterile Box qty
 <p>22F 22M/15F</p>	3504	75	70-150mm	3504000S	75
 <p>22F 22M/15F</p>	3502	75	70-150mm	3502000S	40
 <p>Luer lock port 22F 22M/15F</p>	3514	75	70-150mm	3514000S	40
 <p>Swivel 22F 22M/15F swivel</p>	3512	75	70-150mm		
 <p>'Flip top' cap with 7.6mm port Swivel 22F 22M/15F swivel</p>	3520	75	70-150mm	3520000S	40
 <p>Double 'Flip top' cap with seal Swivel Covers 7.6mm and 9.5mm applications 22F 22M/15F swivel</p>	3521	75	70-150mm	3521000S	40
 <p>Paediatric/neonatal 7.6mm port 15M 15F</p>	3535	15	49-100mm		
 <p>Neonatal 7.6mm port 15M 8.5F</p>	3534	25	49-100mm		