

Thermometer FT 15

Sales Unit: 12

EAN-No.: 42 11125 79420 0

Art.-No.: 794.20



Digital EXPRESS thermometer, extra rapid measurement in only 10 seconds, flexible measuring tip, without glass, mercury free

Type of thermometer:	Digital EXPRESS thermometer
Advantages of use:	Extra rapid measurement in only 10 seconds, flexible measuring tip
Safety:	Without glass, mercury free
Memory:	10 memories
Measurement accuracy:	+/- 0,1 °C
Signal:	Audible signal when measuring stops and if temperature exceeds 38 °C
Large display:	Extra large display
Waterproof:	100% waterproof
Cleaning:	Can be disinfected
Storage case:	With storage case
Battery:	Including battery
Battery changing indicator:	Yes
Automatic switch off:	Yes

755.xxx - 0304 No responsibility accepted for errors