


Bone & Bone Marrow Biopsy

T-Lok™ Bone Marrow Biopsy Needles

T-Lok Bone Marrow Biopsy Needles are designed to provide quality biopsies by entrapping bone marrow in the needle without altering the sample.



Box of 10

Catalogue #	Gauge	Length	Cannula
DBMND0804TL	8G	4" (10.2cm)	w/Extraction Cannula
DBMND0806TL	8G	6" (15.2cm)	w/Extraction Cannula
DBMND1104TL	11G	4" (10.2cm)	w/Extraction Cannula
DBMND1106TL	11G	6" (15.2cm)	w/Extraction Cannula
DBMND1104D	11G	4" (10.2cm)	w/D Shape Extraction Cannula

Bone Marrow Biopsy Needles (Box of 10)

Catalogue #	Gauge	Length
DBMND13025	13G	2.5" (6.4cm)
DBMND1304	13G	4" (10.2cm)

Bone Marrow Aspiration Needles

Box of 10

Catalogue #	Gauge	Length
DBMN1501	15G	2.688" (6.8cm) max.
DBMN1504	15G	4" (10.2cm) max.
DBMN1601	16G	2.688" (6.8cm) max.



Bone Marrow Harvest Needles

The Bone Marrow Harvest Needles are designed for allogenic and autologous bone marrow harvest-transplantation procedures.

Box of 10

Catalogue #	Gauge	Length	Tip
BMHN104VX	11G	4" (10.2cm)	Double Diamond
BMHN104X	11G	4" (10.2cm)	Bevel
BMHN1302X	13G	2" (5.1cm)	Bevel
BMHN1303X	13G	3" (7.6cm)	Bevel
BMHN1304VX	13G	4" (10.2cm)	Double Diamond
BMHN1304X	13G	4" (10.2cm)	Bevel

