

Superset™ catheter mounts

Superset catheter mounts have a concertina design to allow for complete positioning control.

The design allows the adjustment of compressible volume and the adjustment of length from 70mm to 150mm, excluding connectors.

The range is available with fixed and double swivel elbows and a variety of ports for all applications. For added safety all Intersurgical port caps are fixed to the connectors, to prevent loss or misplacement in use.

Please see page 4 for the guide to ports and caps.



Video available at www.intersurgical.com



Code	3504000S
Box Qty.	75
Compressible Volume	22-43ml

Code	3502000S
Box Qty.	40
Compressible Volume	23-44ml

Code	3514000S
Box Qty.	40
Compressible Volume	26-47ml

Code	3520000S
Box Qty.	40
Compressible Volume	23-47ml

Code	3521000S
Box Qty.	40
Compressible Volume	26-50ml

Make an enquiry



Contact us

info@intersurgical.com • www.intersurgical.com