

ARTERIAL CATHETERS



Careflow™ Arterial Catheters

Sterile, single use

Merit Careflow™ Polyurethane catheters are designed for access to radial or femoral arteries, for the determination of blood parameters or monitoring blood pressure.

CAREFLOW™ ARTERIAL CATHETER SETS

PRODUCT CODE	PRODUCT	CATHETER		NEEDLE INTRODUCER		GUIDEWIRE OD & LENGTH	BOX/CASE CONTENT (UNITS)	PRIMING VOLUME (ML)	FLOW RATE (ML/MIN)
		EXTERNAL DIAMETER	LENGTH	GAUGE	LENGTH				
681631	Radial artery set	2.5Fr (22G)	60mm	22G	40mm	0.5 x 250mm	10/100	0.08	12.4
681646	Radial artery set	3Fr (20G)	60mm	22G	40mm	0.5 x 250mm	10/100	0.08	24.0
681642	Femoral artery set	3Fr (20G)	200mm	22G	40mm	0.5 x 450mm	10/100	0.14	15.4
681648	Femoral artery set	4Fr (18G)	200mm	17G	70mm	0.6 x 450mm	10/100	0.24	35.8

Merit Careflow™ Polyurethane catheter set contents:

Merit Careflow™ radiopaque polyurethane catheter, introducer needle, metal guidewire and Luer Lok™ plug.