

ThinPrep® PreservCyt Vials (GYN)



Item No: 70098-002

Item Desc: PRESERVCYT INT'L SOLUTION KIT

Box of 250 vial pack, 20ml/2oz vial

ThinPrep® Filters (GYN)



Item No: 70099-001

Item Desc: FILTER, THINPREP PAP TEST, 500-PACK

Box of 500 filters (5x100 pc)

ThinPrep® Microscope Slides (GYN)



Item No: 70303-001

Item Desc: THINPREP SLIDE KIT

Box of 5 x 100 pc, 500-pack

Slide Kit TP Imaging slide kit



Item No: 70825-001

Item Desc: SLIDE KIT, THINPREP IMAGING, SMS

Box of 5 x 100 pc, 500-pack

CytoLyt® Solution Bottles



Item No: 70408-002

Item Desc: CY TOLYT 32OZ 4-PACK INTERNATIONAL

4 bottle pack (4 x 946ml)
Transport medium and wash solution

PreservCyt® Solution Bottles



Item No: 70406-002

Item Desc: PRESERVCYT 32OZ 4-PACK INTERNATIONAL

4 bottle pack (4 x 946ml)
Fixative solution

Specimen Collection Devices



Item No: 70124-001

Item Desc: ASSEMBLY KIT, MEDSCAND CYTOBRUSH/SPATULA

20 x 25 (500)/box includes: 500 Medscand Cytobrush/Plastic Spatula combination

Specimen Collection Devices



Item No: 70671-001

Item Desc: ROVERS CERVEX-BRUSH KIT

20 x 25 (500)/box

Specimen Collection Devices



Item No: 380101-030

Item Desc: CERVEX BRUSH, GREEN HANDLE, COMBI 20X25 PCS

20 x 25 (500)/box

Dispenser Pump



Item No: 50705-001

Item Desc: DISPENSER PUMP 1-PACK

Use for PreservCyt or Cytolyt
Dispensing 30ml

Stackable ThinPrep® Vial Tray



Item No: 70264-001

Item Desc: SAMPLE VIAL TRAY 4-PACK ASSAY

Box of 4 trays
Each tray can hold 40 vials

ThinPrep® Stains



Item No: 70897-002

Item Desc: THINPREP STAIN SET, INTERNATIONAL

Kit containing 5 stains:

4 liters Rinse Solution (70779-001)
4 liters Nuclear Stain (70780-001)
4 liters Orange G Solution (70781-002)
4 liters EA Solution (70782-002)
4 liters Bluing Solution (70793-001)

ThinPrep® Stains



Item No: 70779-001

Item Desc: THINPREP RINSE SOLUTION,
4L PACKAGED PRODUCT, INTERNATIONAL

ThinPrep® Stains



Item No: 70780-001

Item Desc: THINPREP NUCLEAR STAIN,
4L PACKAGED PRODUCT, INTERNATIONAL

ThinPrep® Stains



Item No: 70781-002

Item Desc: THINPREP ORANGE G SOLUTION
4L PACKAGED PRODUCT, INTERNATIONAL

ThinPrep® Stains



Item No: 70782-002

Item Desc: THINPREP EA SOLUTION 4L
PACKAGED PRODUCT, INTERNATIONAL

ThinPrep® Stains



Item No: 70793-001

Item Desc: THINPREP BLUING SOLUTION
4L PACKAGED PRODUCT, INTERNATIONAL

ThinPrep
PAP TEST

Diagnostic Solutions | Hologic.com | NordicInq@hologic

PB-00747-NOR-EN Rev 003 @2022 Hologic, Inc. All rights reserved. Hologic, The Science of Sure, Cytolyt, ThinPrep, PreservCyt and associated logos are trademarks or registered trademarks of Hologic, Inc. and/or its subsidiaries in the United States and/or other countries. All other trademarks, registered trademarks and product names are the property of their respective owners. This information is intended for medical professionals and is not intended as a product solicitation or promotion where such activities are prohibited. Because Hologic materials are distributed through websites, eBroadcasts and tradeshows, it is not always possible to control where such materials appear. For information on specific products available for sale in a particular country, please contact your Hologic representative or write to NordicInq@hologic

The content of this text is for information purposes and is not intended to be medical advice. Contact your doctor for specific advice on your health and treatment.