



## Micro pipettes, end-to-end, Vitrex® Micro

Leverantör: [Vitrex Medical A/S](#)



For single use.

- Volume is determined by both ends
- Packed in glass vials
- Vials have a dispenser top, making it easy to shake out one tube at a time, therefore ensuring hygienic handling
- Dispenser tops are colour coded to indicate the tube volume in conformity with ISO 1769

Accuracy:  $\pm 0,50\%$  (0,5 to 10  $\mu\text{l}$ );  $\pm 0,50\%$  (20 to 100  $\mu\text{l}$ )  
 CV: 0,75% (0,5 to 10  $\mu\text{l}$ ); 0,60% (20 to 100  $\mu\text{l}$ ).

**Förpackning:** Ten vials are packaged in a sturdy carton in order to protect the tubes from damage during handling or transport.

Volym	Längd	Typ	Förp.	VWR art.nr
10 $\mu\text{l}$	40 mm	Plain	1.000	MODU165513
20 $\mu\text{l}$	40 mm	Plain	1.000	MODU165613
25 $\mu\text{l}$	50 mm	Plain	1.000	MODU165713
50 $\mu\text{l}$	67 mm	Plain	1.000	MODU165913
100 $\mu\text{l}$	100 mm	Plain	1.000	MODU166013

