



Product description

Order number	95.1123
Product description	Mailing container for cool transport, english

Product characteristics

Label/ Print	With printed label
Colour of print/label	Red
Design	English cool transport container
Example of use	Cool transport
Absorbent liner	no

Size

Width of product	71 mm
Height of product	170 mm
Length of product	140 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Acrylnitrile-butadiene-styrole (ABS)
Colour of product	White

Purity & certification

Satisfies the requirement	P 650, ADR
Product in accordance with MDD or IVDD	No marking necessary
CE certified	No CE marking
Lot or serial number	batched

Packaging

Minimum order qty.	6
Pcs. / inner box	1
Pcs. / box	6
Pcs. / pallet	648
Depth of case	308 mm
Width of case	223 mm
Height of case	195 mm
Case volume	0,0134 cbm
Weight of case	2,3 KG
EAN of inner box	4038917085403
EAN of outer box	4038917004084

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

