



### **Product description**

|                     |                                                                                                                                        |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| Order number        | 57.462                                                                                                                                 |
| Product description | Tube, working volume: 12 ml, (LxØ): 110 x 17 mm, material: PS, conical base, transparent, push cap, with graduations, 500 piece(s)/bag |

### **Product characteristics**

|                          |              |
|--------------------------|--------------|
| Graduation               | yes          |
| Closure type             | push cap     |
| Base shape               | conical base |
| Volume of work           | 12 ml        |
| Centrifugation max (RCF) | 2500 x g     |

### **Size**

|                   |        |
|-------------------|--------|
| Diameter          | 17 mm  |
| Length of product | 110 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

|                   |                  |
|-------------------|------------------|
| Product material  | Polystyrene (PS) |
| Colour of product | transparent      |

### **Purity & certification**

|                           |                                |
|---------------------------|--------------------------------|
| Satisfies the requirement | sealed at 0.6 bar/10 min       |
| Product category          | In vitro diagnostic device, CE |
| CE certified              | CE                             |
| Batched                   | yes                            |

### **Packaging**

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 1000          |
| Type of smallest subpackaging | bag           |
| Piece(s) / inner box          | 500           |
| Piece (s) / outer case        | 1000          |
| Piece(s) / pallet             | 36000         |
| Depth of case                 | 398 mm        |
| Width of case                 | 398 mm        |
| Height of case                | 300 mm        |
| Case volume                   | 0.0475 cbm    |
| Weight of product             | 0.0063 kg     |
| Weight of case                | 6.98 kg       |
| EAN of inner box              | 4038917081702 |
| EAN case                      | 4038917000819 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

