



Product description

Order number	76.1414.100
Product description	Lid without spout for urine container 125 ml, 100 pieces packed

Product characteristics

Type of cap	Cap
Design	Without spout
Application	Collection and storage of urine
Example of use	Cap for urine collection tube
Suitable for	75.1355 75.1355.100

Material & colours

Product material	Polyethylene (PE)
Colour of product	Red

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product in accordance with MDD or IVDD	IVD
CE certified	CE
Lot or serial number	batched

Packaging

Minimum order qty.	1000
Pcs. / inner box	100
Pcs. / box	1000
Pcs. / pallet	24000
Depth of case	590 mm
Width of case	390 mm
Height of case	265 mm
Case volume	0,061 cbm
Weight of case	3,84 KG
EAN of inner box	4038917128988
EAN of outer box	4038917417068

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

