



### PVC Catheter Mounts with 90° Fixed Elbow

Description	Qty/box	REF
10 cm, 22 F adapter -22M/15 F patient connector, sterile	25	330/5109
15 cm, 22 F adapter -22M/15 F patient connector, sterile	25	330/5110
10 cm, 15M adapter -22M/15 F patient connector, sterile	25	330/5307
15 cm, 15M adapter -22M/15 F patient connector, sterile	25	330/5308



### PVC Catheter Mounts with 90° Elbow with CO<sub>2</sub> Sampling Port

Description	Qty/box	REF
10 cm, 15M adapter -22M/15 F patient connector, sterile	25	330/5321
10 cm, 15M adapter -22M/15 F patient connector with cap tethered on port, sterile	25	330/5321TC
15 cm, 15M adapter -22M/15 F patient connector, sterile	25	330/5322
15 cm, 15M adapter -22M/15 F patient connector with cap tethered on port, sterile	25	330/5322TC
10 cm, 22 F adapter -22M/15 F patient connector, sterile	25	330/5645
15 cm, 22 F adapter -22M/15 F patient connector, sterile	25	330/5646
15 cm, 22 F adapter -22M/15 F patient connector with cap tethered on port, sterile	25	330/5646TC



### PVC Straight Catheter Mounts

Description	Qty/box	REF
10 cm, 15M-15 F, sterile	25	332/5114
15 cm, 15M-15 F, sterile	25	332/5115
10 cm, 22 F-15 F, sterile	25	332/5309
15 cm, 22 F-15 F, sterile	25	332/5310
5 cm, 22 F-15 F, sterile	25	332/5316
5 cm, 15M-15 F, sterile	25	332/5317
20 cm, 15M-15 F, sterile	25	332/5332
20 cm, 15M adapter – 22M/15 F patient connector, sterile	25	332/5685
30 cm, 15M adapter – 22M/15 F patient connector, sterile	25	332/5339



### Extendible Catheter Mounts with Double Swivel Elbow With Suction and Bronchoscopy Caps

Description	Qty/box	REF
7/16 cm, 15M adapter – 22M/15 F patient connector, sterile	25	331/5661
7/16 cm, 22 F adapter – 22M/15 F patient connector, sterile	25	331/5669

