



Product Name	Oxygen Thick Wall Tube	
Reference numbers	1180000	
Product Type	Oxygen Tubing	
Manufacturer	Intersurgical	
Materials		PVC
Packaging		
Unit Packaging	Individually packaged	
Box quantity	1	
Pallet quantity	72	
Instructions for Use	Instructions for use with product	No
	Instructions for use with each box	No
CE Classification	Class 2A	
Applicable Standards	BS EN ISO 13485: Medical devices — Quality management systems — Requirements for regulatory purposes	
	ISO 9001: Quality Management Systems - Requirements	
	ISO 15223-1: Medical Devices - Symbols to be Used with Medical Device Labels, Labelling, & Information to be Supplied	
	BS EN 13544-2: Respiratory therapy equipment. Tubing and connectors	
Sterile	No	
Shelf life	5 Years	
Disposal Methods	As hospital guidelines	
Re-use		Single use only

Product characteristics

Technical information

Roll Length (m)	50
Maximum outside diameter (mm)	8
Minimum inside diameter (mm)	5



Material Abbreviations

ABS	Acrylonitrile Butadiene Styrene
AL	Aluminium
CR	Polychloroprene Rubber
EDPM	Ethylene Propylene Terpolymer Rubber
HDPE	High Density Polyethylene
HIPS	High Impact Polystyrene
LDPE	Low Density Polyethylene
LSR	Liquid Silicone Rubber
MBS	Methacrylate Butadiene Styrene
NRL	Natural Rubber Latex
PA	Polyamide (Nylon)
PC	Polycarbonate
PE	Polyethylene
POM	Polyoxymethylene
PP	Polypropylene
PS	Polystyrene
PSU	Polysulfone
PU	Polyurethane
PVC	Poly Vinyl Chloride
SB	Styrene Butadiene
Si	Silicone
SS	Stainless Steel
TPE	Thermoplastic Elastomer
SBC	Styrene Butadiene Copolymer

