



Careflow™ Central Venous Catheters

Single Lumen without Extension

Sterile, single use

- Careflow™ is a range of paediatric and adult polyurethane central venous catheters designed for catheterisation via the Seldinger technique.
- Range of sets contains a Nitinol guidewire and the Venaguide™ Seldinger wire-advancing device facilitating single-handed guidewire insertion.
- Tapered atraumatic catheter tip designed to reduce the risk of venous damage while improving catheter insertion.
- The streamlined junction hub enables ease of catheter suture and comfortable positioning on the patient's skin.
- Careflow™ catheters offer 1cm and 5cm interval depth markings to facilitate catheter placement and clear lumen size and hub colour coding indication.
- Floswitch® and Floswitch® Introducer Catheter options available throughout range.
- Pressure rated with a maximum pressure tolerance of 45 psi.

CAREFLOW™ SINGLE LUMEN SETS WITHOUT EXTENSION

CODE WITHOUT FLOSWITCH®	CODE WITH FLOSWITCH®	EXTERNAL DIAMETER	LENGTH	CANNULA INTRODUCER		FLOSWITCH® INTRODUCER		NEEDLE INTRODUCER		GUIDEWIRE OD & LENGTH	SECONDARY FIXATION (SFD)	BOX/CASE CONTENT (UNITS)	FLOW RATE (ML/MIN)
				GAUGE	LENGTH	GAUGE	LENGTH	GAUGE	LENGTH				
681614	681623	2.5Fr (22G)	100mm	N/A	N/A	20G	45mm	22G	40mm	0.5 x 350mm		10/100	10
681639	681640 [≠]	3Fr (20G)	150mm	N/A	N/A	20G	45mm	22G	40mm	0.5 x 350mm		10/100	16
681643	681644	3Fr (20G)	200mm	N/A	N/A	20G	45mm	22G	40mm	0.5 x 450mm		10/100	15
681649	681650	4Fr (18G)	200mm	N/A	N/A	18G	60mm	17G	70mm	0.6 x 450mm		10/100	36
681669	681670	5Fr (16G)	150mm	N/A	N/A	15G	70mm	17G	70mm	0.9 x 450mm	Yes	10/100	71
681672	681673	5Fr (16G)	200mm	N/A	N/A	15G	70mm	17G	70mm	0.9 x 450mm	Yes	10/100	63
681706	681707	7Fr (14G)	150mm	N/A	N/A	14G	70mm	17G	70mm	0.9 x 450mm	Yes	10/100	190
681709	681710	7Fr (14G)	200mm	N/A	N/A	14G	70mm	17G	70mm	0.9 x 450mm	Yes	10/100	136

PRODUCT CONTENTS

Single lumen Careflow™ with Floswitch® set contents:

Careflow™ Polyurethane catheter, steel needle introducer, Floswitch® overneedle introducer, vessel dilator, Venaguide™ Seldinger wire-advancing device with Nitinol guidewire, Floswitch® revolving Luer-Lok™ system, luer lock plug and secondary fixation device (where shown).

Single lumen Careflow™ without Floswitch® set contents:

Careflow™ Polyurethane catheter, steel needle introducer, vessel dilator, Venaguide™ Seldinger wire-advancing device with Nitinol guidewire and secondary fixation device (where shown).

[≠] Note: Part number 681640 contains 0.5 x 450mm guidewire.