


PRODUKTINFORMATION-TRAKEALTUB



IVORY PVC, ORAL/NASAL

- Knick-resistent
- Utan kuff
- Tydliga markeringar
- 15 mm ISO standard konnektor
- Mjuka avfasade kanter vid tubspets



Benämning	Storlek (mm)	Vårt Art nr	Avdfp/st
Trakealtub, IVORY utan kuff	2,5	100/105/025	10
Trakealtub, IVORY utan kuff	3,0	100/105/030	10
Trakealtub, IVORY utan kuff	3,5	100/105/035	10
Trakealtub, IVORY utan kuff	4,0	100/105/040	10
Trakealtub, IVORY utan kuff	4,5	100/105/045	10
Trakealtub, IVORY utan kuff	5,0	100/105/050	10
Trakealtub, IVORY utan kuff	5,5	100/105/055	10
Trakealtub, IVORY utan kuff	6,0	100/105/060	10
Trakealtub, IVORY utan kuff	6,5	100/105/065	10
Trakealtub, IVORY utan kuff	7,0	100/105/070	10
Trakealtub, IVORY utan kuff	7,5	100/105/075	10
Trakealtub, IVORY utan kuff	8,0	100/105/080	10
Trakealtub, IVORY utan kuff	8,5	100/105/085	10
Trakealtub, IVORY utan kuff	9,0	100/105/090	10