

KLINIFIX - COHESIVE FIXATION BANDAGE

Description

Cohesive elastic securing bandage, white. Latex free.

Product information

In compliance with: MDR (EU) 2017/745 Class I



Article overview

Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
3016838	KLINION KLINIFIX COHESIVE ELASTIC FIXATION BANDAGE SELF-ADHESIVE LATEX FREE 4 M X 2,5 CM REF 132377	10 pcs.	450 pcs.	8715343029671	8715343029688	132377
3016830	KLINIFIX COHESIVE ELASTIC BANDAGE LATEX FREE 4MX4CM 132378	pcs.	pcs.			132378
3038996	KLINION KLINIFIX COHESIVE ELASTIC FIXATION BANDAGE LATEX FREE 4 M X 6 CM REF 132379	pcs.	pcs.			132379
3038997	KLINION KLINIFIX COHESIVE ELASTIC FIXATION BANDAGE LATEX FREE 4 M X 8 CM REF 132380	pcs.	pcs.			132380
3039000	KLINION KLINIFIX COHESIVE ELASTIC FIXATION BANDAGE LATEX FREE 4 M X 10 CM REF 132380	pcs.	pcs.			132381
3039001	KLINIFIX COHESIVE ELASTIC BANDAGE LATEX FREE 4MX12CM 132382	5 pcs.	140 pcs.	8715343060292	8715343051665	132382
3016835	KLINIFIX COHESIVE ELASTIC BANDAGE LATEX FREE 20MX6CM 132386	1 pcs.	24 pcs.	8715343060346	8715343008355	132386
3016836	KLINION KLINIFIX COHESIVE ELASTIC FIXATION BANDAGE SELF-ADHESIVE LATEX FREE 20 M X 8 CM REF 132387	pcs.	pcs.			132387
3016837	KLINION KLINIFIX COHESIVE ELASTIC FIXATION BANDAGE SELF-ADHESIVE LATEX FREE 20 M X 10 CM REF 132388	pcs.	pcs.			132388