

White colour
Cor branca



Tubular bandage

Ligadura tubular

Bandage tubulaire


Schlauchverband

Buisvormig verband

Benda tubolare

Venda tubular

Tubförfband



TUBULAR BANDAGES SEAMLESS

- Viscose.
- Protection of several parts of the body.
- Ideal for use next to the skin.
- Seamless and smooth. Easy application even on body parts that are difficult to dress. Run off proof in any direction. Sterilizable.
- Stretchable widthwise.
- Stretched length: 10m and 20m.
- Dispenser shelf box + transport carton.

LIGADURAS TUBULARES, SEM COSTURA

- Viscose.
- Proteção de várias partes do corpo.
- Ideal para uso junto à pele.
- Sem costuras e macias. Fácil aplicação mesmo nas partes mais difíceis do corpo. Indesmalhável em qualquer direcção. Esterilizáveis.
- Extensíveis à largura.
- Comprimento esticado: 10m e 20m.
- Caixa de cartolina dispenser + caixa de transporte.

Non sterile • Não esterilizada

REF	Size Tamanho	Pcs./Sh.box Uni./Cx.cartolina	Pcs./Tr.carton Uni./Cx.transporte
4211-001	20 m x 1,5 cm	1	20
4211-002	20 m x 2,5 cm	1	12
4211-003	20 m x 3,5 cm	1	10
4211-004	20 m x 5,5 cm	1	6
4211-005	20 m x 6,5 cm	1	5
4211-006	20 m x 7 cm	1	5
4211-007	20 m x 9 cm	1	4
4211-008	20 m x 12 cm	1	3
4211-009	10 m x 16 cm	1	18
4211-010	10 m x 21 cm	1	8
4211-011	20 m x 4 cm	1	8
4211-012	20 m x 6 cm	1	6
4211-013	20 m x 8 cm	1	4
4211-014	20 m x 10 cm	1	3