

## Anti-Siphon valve

04302260452



- Fyllnadsvolym 0.2 ml
- 1 Hon-luer
- 1 Anti-sifonventil
- 1 Han-luer
- Lipidresistent
- Latexfri
- PVC-fri
- Kvantitet 300/100x3



300/100x3



1 Han-luer



1 Hon-luer



1 Anti-sifonventil



Latexfri



DEHP-fri



PVC-fri



0086

Customer Service  
+ 46 (0) 8 5444 3200  
(0) 8 5444 3225  
[CS.Sweden@carefusion.com](mailto:CS.Sweden@carefusion.com)

Carefusion Sweden 314 AB  
Sollentuna