

ELECTROSURGERY

Electrosurgical Electrodes, Sterile, Single Use Stainless Steel Electrodes



E1550 150 UNITS/CASE

Stainless Steel Ball Electrode

- Total length: 5,33 cm (2,1 in.)
- Active length: 1 cm (0,4 in.)
- Diameter: 5 mm (3/16 in.)



E1552 150 UNITS/CASE

Stainless Steel Needle Electrode

- Total length: 7,2 cm (2,84 in.)
- Active length: 2,8 cm (1,1 in.)
- Tip radius: 0,762 mm (0,030 in.)



E1551G 150 UNITS/CASE

Stainless Steel Blade Electrode

- Total length: 6,2 cm (2,44 in.)
- Active length: 2,8 cm (1,1 in.)



E1552-6 50 UNITS/CASE

Stainless Steel Extended Needle Electrode

- Total length: 16,51 cm (6,5 in.)
- Active length: 2,54 cm (1 in.)
- Tip radius: 0,762 mm (0,030 in.)



E1551X 150 UNITS/CASE

Stainless Steel Blade Electrode

- Hex locking
- Total length: 6,2 cm (2,44 in.)
- Active length: 2,54 cm (1 in.)



E2401 100 UNITS/CASE

LectroBrasive™ Tip Cleaner

For use with stainless steel electrodes

- Sterile, Single Use
- 5,08 cm × 5,08 cm (2 in. x 2 in.)
- Adhesive backing
- X-ray detectable



E1551-6 50 UNITS/CASE

Stainless Steel Extended Blade Electrode

- Total length: 16,51 cm (6,5 in.)
- Active length: 2,54 cm (1 in.)