



Product description

Order number	95.64.250
Product description	Styrofoam container, white, for tubes Ø 15 and 16 mm, can be parted in the middle at perforation, External dimensions (WxLxH): 107 x 416 x 50 mm, grid: 5 x 20, for 100 tubes

Product characteristics

Suitable for	Tubes Ø 15 and 16 mm
Storage capacity	100
Format	5 x 20

Size

Width of product	107 mm
Height of product	50 mm
Length of product	416 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material	Styrofoam
Colour of product	White
Cap material	Styrofoam
Colour of cap	White

Purity & certification

Product in accordance with MDD or IVDD	No marking necessary
CE certified	No CE marking
Lot or serial number	non batched

Packaging

Minimum order qty.	1
Pcs. / inner box	1
Pcs. / box	28
Pcs. / pallet	448
Depth of case	498 mm
Width of case	428 mm
Height of case	380 mm
Case volume	0,0809 cbm
Weight of case	1,65 KG
EAN of inner box	4038917095587
EAN of outer box	4038917007764

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

