



Product description

| | |
|---------------------|--|
| Order number | 95.902 |
| Product description | Large post-transport packaging for diagnostic samples in accordance with packaging regulation P 650, oversize letter, structural authorisation, flat design (length 220 mm x width 170 mm x height 40 mm), 50 sterile pieces packed non-erected in the bag |

Product characteristics

| | |
|----------------|--|
| Label/ Print | With printed label |
| Example of use | External packaging according to P 650/ ADR |

Size

| | |
|-------------------|--------|
| Width of product | 170 mm |
| Height of product | 45 mm |
| Length of product | 220 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Material & colours

Product material Cardboard

Purity & certification

Satisfies the requirement P 650, IATA, ADR
Product in accordance with No marking necessary
MDD or IVDD
CE certified No CE marking
Lot or serial number batched

Packaging

Minimum order qty. 50
Pcs. / inner box 50
Pcs. / box 50
Pcs. / pallet 2400
Depth of case 488 mm
Width of case 401 mm
Height of case 145 mm
Case volume 0,0284 cbm
Weight of case 3,86 KG
EAN of inner box 4038917084574
EAN of outer box 4038917069854

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

