





## Alaris® GP Series Primary infusion set with 15µm chamber filter

60593



- 265 cm Aggregatlängd
- Fyllnadsvolym 24 ml
- Bore 3 mm Diameter
- 1 Droppkammare med spike (ventilerad (e)/oventilerad(e))
- Droppar per ml 20
- Hydrofobiskt luftfilter (med lock)
- 1 Rullklämma
- Alaris® säkerhetsklämma
- Pumpsegment av silikon
- 1 Han-luer
- Polyvinylklorid (PVC)
- Latexfri
- Kvantitet 80



15 µm Droppkammfilter



1 Han-luer



Latexfri



DEHP-fri

Carefusion Sweden 314 AB  
Sollentuna