



Tejp Sterilt grön/svart 25mmx66m

General data

Customer ID : 100001
Contact :
Contact email :
Representative :
Job Manager : Admin
Internal contact : Admin
Ref. at customer : MT115
Press :

Status data

Prepress : Approval OK
Status : Utan ändring
Korrektur : Inget korr

Technical data

Size : 25,000 mm (W) x 66,000 mm (H)
Shape :
Tool : N/A
Substrate : N/A
Description : N/A
Front printing : N/A
Adhesive printing : N/A
Back printing : N/A

Colours

Front  Pantone Green -----
Front  Black -----

Texts / Comments / Info

Product info :
Region Östergötland Fraktfritt
Region Östergötland Fraktfritt
Region Östergötland Fraktfritt
Prepress instructions : N/A
Remarks : N/A

Packing