



### Product description

Order number	75.562.105
Product description	Container with screw cap 100 ml, with graduation, PP, transparent, with assembled lid, length 73 mm, opening diameter 62 mm, sterile, with safety label, 5 pieces in the bag

### Product characteristics

Label/ Print	With safety label
Graduation	Yes
Type of cap	Screw cap
Design	Container with assembled lid, sterile, with safety label
Application	Collection and storage of urine
Example of use	Collection tube for urine samples
Suitable for	76.564

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Size**

Sample volume	100 ml
Max. volume	100 ml
Diameter	57 mm
Diameter of opening	62 mm
Height of product	76 mm
Length excluding cap	72 mm

### **Material & colours**

Product material	Polypropylene (PP)
Colour of product	Transparent
Cap material	High Density Polyethylene (HD-PE)

### **Purity & certification**

Product in accordance with MDD or IVDD	IVD
CE certified	CE
Sterilisation	sterile
Lot or serial number	batched

### **Packaging**

Minimum order qty.	200
Type of packaging	5 sterile pieces packed
Pcs. / inner box	5
Pcs. / box	200
Pcs. / pallet	5600
Depth of case	588 mm
Width of case	398 mm
Height of case	255 mm
Case volume	0,0597 cbm
Weight of case	3,78 KG
EAN of inner box	4038917113632
EAN of outer box	4038917014052

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

