

# Produktblad

## Tallrik djup 480 ml

### BESKRIVNING

Tallrikar i fibermaterialet bagasse - en restprodukt från sockerproduktionen.



### Produktinformation

Artikelnummer	3714114
Beskrivning	Tallrik djup 480 ml
Förpackningsformat/mått	Ø 17,5x4 cm
Färg	Naturvit
Försäljningsenhet	TFP
Vikt försäljningsenhet	12,00 kg
Ean säljenhet	
Förpackningsavgift	Ja
Punktskatt	Nej
Antal/Innerförpackning	50 st
Antal/TFP	1000 st

### Materialinformation

Papper/Kartong/Bagasse/Pulp

[Materialguide med de vanligast förekommande materialen i våra take away-produkter. \(mediapaper.se\)](#)

 Artikeln tål uppvärmning i mikrovågsugn.

### GODKÄND FÖR LIVSMEDELSKONTAKT

We certify that these products fulfill the requirements on products intended for use in contact with food as describe in following European legislation:

- EU Framework Regulation EC No 1935/2004 2004-10-27, relating to requirement for all food contact material amending Directive 80/590/ EC, 89/109/EC
- Regulation EC No 2023/2006 2006-12-22, relating to good manufacturing practice for food contact materials with later modification
- LFGB (BfrXXXVI and BfrXXXVI/2)



PAC-PRODUCTION SWEDEN AB  
BERGLUNDAVÄGEN 5-7  
BOX 409  
701 48 ÖREBRO  
kontakt@pac.se

Telefon  
019-175700  
Faxnr  
019-175753