


Product Name	22mm Flextube - 50m roll cuffed at 400mm	
Reference numbers	1574000	
Product Type	Flextube™	
Manufacturer	Intersurgical	
Materials	LDPE	
Packaging	Individually packaged	
Unit Packaging	1	
Box quantity	12	
Pallet quantity		
Instructions for Use	Instructions for use with product	N/A
	Instructions for use with each box	N/A
Certification	CE0120	Class 2A
Applicable Standards	BS EN ISO 13485: Medical devices — Quality management systems — Requirements for regulatory purposes	
	ISO 9001: Quality Management Systems - Requirements	
	ISO 15223-1: Medical Devices - Symbols to be Used with Medical Device Labels, Labelling, & Information to be Supplied	
	BS EN ISO 5367: Anaesthetic and Respiratory Equipment. Breathing Sets and Connectors	
Sterile	No	
Shelf life	5 years	
Disposal Methods 	As hospital guidelines	
Re-use	Single use only	

Product characteristics

Roll length (m)	50
Cuff length (mm)	400

Technical information - Results per m or 22mm flextube

Resistance to flow (mbar) at 60l/min	0.2
Compliance at 100cmH2O (ml/cmH2O/m)	0.5
Compressible volume (ml/m)	504



**Material
Abbreviations**

ABS	Acrylonitrile Butadiene Styrene
AL	Aluminium
HDPE	High Density Polyethylene
LDPE	Low Density Polyethylene
PC	Polycarbonate
PE	Polyethylene
PP	Polypropylene
PVC	Poly Vinyl Chloride
SB	Styrene Butadiene
Si	Silicone
SS	Stainless Steel
TPE	Thermoplastic Elastomer

