



### **Product description**

Order number	95.1706
Product description	Wall holder made of stainless steel, for Formalin System, 1 pc./case

### **Product characteristics**

Application	Capillary and venous blood collection
-------------	---------------------------------------

### **Size**

Width of product	101 mm
Length of product	200 mm

### **Material & colours**

Product material	Stainless steel
------------------	-----------------

### **Purity & certification**

Product in accordance with MDD or IVDD	IVD
CE certified	CE
Lot or serial number	batched

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





### **Packaging**

Minimum order qty.	1
Type of packaging	Case
Pcs. / inner box	1
Pcs. / box	1
Pcs. / pallet	324
Depth of case	268 mm
Width of case	191 mm
Height of case	100 mm
Case volume	0,005 cbm
Weight of case	0,42 KG
EAN of inner box	4038917335461
EAN of outer box	4038917335454

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

