



Product description

Order number	09.1704.001
Product description	S-Monovette® 9 ml, for Formalin System, 92 x 16 mm, Paper label, 50 pcs./bag

Product characteristics

Type of preparation	Without preparation
Label/ Print	With paper label
Application	Capillary and venous blood collection

Size

Diameter	16 mm
Length of product	92 mm

Material & colours

Product material	Polypropylen (PP)
Membrane material	Cautchouc
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





Purity & certification

Product in accordance with MDD or IVDD	IVD
CE certified	CE
Lot or serial number	batched

Packaging

Minimum order qty.	500
Type of packaging	Bag
Pcs. / inner box	50
Pcs. / box	500
Pcs. / pallet	24000
Depth of case	388 mm
Width of case	378 mm
Height of case	215 mm
Case volume	0,0315 cbm
Weight of case	4,4 KG
EAN of inner box	4038917335317
EAN of outer box	4038917335300

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

