

## Wheelchair Scale MBWA-300 (Class III MDD)

This wheelchair scale with robust handrail allows weighing of wheelchair bound patients, patients standing with a walker and patients sitting on a stationary chair. Durable, large 100 x 94 cm platform, and adjustable rubber feet for stability. The sturdy platform has a low angle, and integrated ramps for easy-on, easy-off access.

The scale can be quickly folded to reduce storage space. Built-in handle and transport wheels for easy movement.

### Technical Data

- Capacity 300 kg
- Division 100 g
- Total Weight 33,8/ kg
- Product Size 1250 x 940 x 1150 mm
- Platform Size 1000 x 940 mm (flat)

### Features

- Movement compensation for easy reading
- BMI function
- Zero/Hold/Print/kg/lb
- 25 mm large LCD with backlight
- Adjustable display angle
- AC adaptor or AA batteries (not included)



**“INNOVATIVE APPLICATIONS  
LATEST TECHNOLOGY  
LONG EXPERIENCE”**

Specifications are subject to change without notice