


"Quickstart" Universal Polystyrene Containers 30ml

Product Code	128A*	128B*	128C*	128BBAC*		
Description	No Label	Printed Label	Plain Label	Boric Acid, Printed Label		
1	Dimensions	OD (mm)	Height (mm)			
		31	94			
2	Packaging	Bag Quantity		Carton Quantity		
		50		400		
3	Material composition	Container: General Purpose Polystyrene (PS) Cap: High density polyethylene (HDPE)				
4	Sterility	Non-sterile Manufactured in an ISO Class 7 environment Clean Room				
5	Quality	Leak tested in accordance with BS EN ISO 6717:2021 Suitable for centrifugation at 3,800 x g 95kPa compliant at ambient temperatures Contains boric acid (128BBAC) - 0.4g ±0.02g				
6	Autoclavability	Not suitable				
7	Chemical Resistance	Dilute acids, bases, alcohols, oils & minerals				
8	Shelf Life	5 years from date of manufacture				
9	CE Marking	Not Applicable				