



GKE-översikt



Lastkontroll	Art.nr.	Bowie-Dick test	Art.nr.
Orange	10011841	Blå	10012147
Startset (inkl. 100 indikat.)	10011975	Startset (inkl. 100 indikat.)	10012146
Grön	10010762	Lila	10012145
Startset (inkl. 100 indikat.)	10011199	Startset (inkl. 100 indikat.)	10011113
Indikatorer 500 st	10011974	Indikatorer 250 st	10011847
Plastlock svart, M10 (mm Ø)	10011845	Plastlock svart, M10 (mm Ø)	10011845
Plastlock svart, M12 (mm Ø)	10011844	Plastlock svart, M12 (mm Ø)	10011844
O-ring / packning, svart 1 st M10 + 1 st M12 / förp.	10011843	O-ring / packning, svart 1 st M10 + 1 st M12 / förp.	10011843
Indikatorhållare	10011842	Indikatorhållare	10011848



Dokumentationssystem	Art.nr.	Kontaineretikett 7 x 2,9 cm	Art.nr.
Märkpistol	10017571	Tapplängd 7 cm	10011976
Etikett med indikator	10010745/ 10013655 /10013657		
Etikett utan indikator	10010743/10010744 / 10013654		
Färgdyna	10010746		