

Product Specifications

Brand : Klinion® Protection
Product : Disposable Aprons

Packaging & Product Reference

Product REF	Description	Length in CM	Width in CM	Colour	Thickness in MM	Packaging	Pieces per Inner carton	Pieces per Outer Carton
522700	Klinion Protection Apron	117	71	White	0,018	Folded in dispenser	100	500
522700-1	Klinion Protection Apron	117	71	White	0,023	Folded in dispenser	100	500
522701	Klinion Protection Apron	125	81	Transparent	0,0195	Folded in dispenser	50	500
522701-1	Klinion Protection Apron	125	81	Transparent	0,025	Folded in dispenser	50	500
522702	Klinion Protection Apron	125	81	Transparent	0,0195	On roll in dispenser	50	500
522703	Klinion Protection Apron	160	90	White	0,019	Folded in dispenser	70	420
522704	Klinion Protection Apron	140	80	White	0,031	Folded in dispenser	70	420
522705	Klinion Protection Short Sleeve Apron	120	98	White	0,03	Folded. Single packed in dispenser	30	300
522706	Klinion Protection Long Sleeve Apron	120	98	Blue	0,03	Folded. Single packed in dispenser	20	100
522706-1	Klinion Protection Long Sleeve Apron	140	98	Blue	0,03	Folded. Single packed in dispenser	15	75
522707	Klinion Protection Apron	150	80	Transparent	0,0195	On roll in dispenser	50	500
522708	Klinion Protection Apron	150	80	Transparent	0,0195	Folded in dispenser	40	400
522709	Klinion Protection Apron	160	90	White	0.03	Folded in dispenser	70	420
522110	Klinion Protection Apron	110	80	White	0.018	Folded in dispenser	100	1000
522130	Klinion Protection Apron	125	80	White	0.018	Folded in dispenser	100	1000
522140	Klinion Protection Apron	140	80	White	0.018	Folded in dispenser	80	800
522160	Klinion Protection Apron	165	80	White	0.018	Folded in dispenser	70	700

* Sizes are approximate

Product Specifications

Product description


Disposable aprons, made of Polyethylene and available in different sizes, colours & dimensions

Intended use

Klinion Personal Protection aprons are intended to be used as barrier protection against harmless fluids and dirt. Protection against general and minimal risks only, which do not expose user to temperatures > 50°C or to dangerous impact

User Info available at www.medeco.org

Dimensions & Design



Features

- Easy fit with two back ties and head loop
- Single-use and disposable
- Available in various sizes to provide sufficient length and width to permit comfortable flexing to accommodate different applications

Product components

- Polyethylene
- Products listed on this specification are:
 - Latex-free
 - DEHP-free
 - Absent of human or animal tissue
 - Absent of blood or blood derivative

Product Specifications

Classification

Category I under Personal Protective Equipment Regulation (EU) 2016/425

Sterilization

N/A – Product is provided non-sterile

Storage conditions

Store at room temperature. Opened boxes should be kept away from fluorescent and direct sunlight. Aprons are packed in dispenser which is suitable for transport. Keep the aprons in the box when not in use.

Shelf life


522700-522709: 3 years

522110-522160: 5 years

Pictures



Ref. product 522700




Ref. product 522707



Ref. product 522705



Ref product 522706

Approval	Name	Date	Signature
Product Manager	Levine van Buren	17-1-2019	
RA-Manager	Francis Dekker	17-1-2019	