

## KLINITULLE - Ointment Compress

### Description

Paraffin compress. Nonwoven impregnated with white paraffin. For superficial wounds and burns. Change daily. Cover with secondary bandage. Medical device. Sterile if pouch and seal are intact.

### Product information

In compliance with: MDD (EU) 93/42/EEC Class IIa



## Article overview

Article number	Description	Inner carton	Outer carton	EAN (IC)	EAN (OC)	Reference
3038895	KLINION KLINITULLE OINTMENT DRESSING 5X5CM 170097 STERILE	10 pcs.	240 pcs.			170097
3035382	KLINION KLINITULLE OINTMENT DRESSING 10 X 10 CM PACKED PER 10 PCS REF 170098 *S*	10 pcs.	240 pcs.			170098
3016980	KLINION KLINITULLE OINTMENT DRESSING 10 X 30 CM REF 170099 *S*	20 pcs.	200 pcs.			170099