

## KLINION - APRONS

### Description

Disposable apron made of polyethylene (PE) material. Barrier protection against harmless fluids and dirt. Protection against general and minimal risks only, which do not expose user to temperatures > 50°C or to dangerous impacts. Before usage inspect the apron for any defect or imperfections. For single use only. Re-use can cause contamination.


### Product information

In compliance with: PPE Regulation (EU) 2016/425 Category I

- Available in many sizes and thickness
- Also available in short and long sleeve



### Article overview

|   | Article number | Description  | Inner carton | Outer carton | EAN (IC)      | EAN (OC)      | Reference |
|---|----------------|--|--------------|--------------|---------------|---------------|-----------|
|  | 3025536        | KLINION PROTECTION APRON PE 110 X 80 CM 100 PCS WHITE REF 522110                   | 100 pcs.     | 1000 pcs.    | 8715343034217 | 8715343034231 | 522110    |
|  | 3025537        | KLINION PROTECTION APRON PE 125 X 80 CM 100 PCS WHITE REF 522130                   | 100 pcs.     | 1000 pcs.    | 8715343034316 | 8715343034330 | 522130    |
|  | 3025538        | KLINION PROTECTION APRON PE WHITE 140CM BOX PE 140 X 80 CM 80 PCS WHITE REF 522140 | 80 pcs.      | 800 pcs.     | 8715343034378 | 8715343034385 | 522140    |
|  | 3025539        | KLINION PROTECTION APRON PE 165 X 80 WHITE REF 522160                              | 70 pcs.      | 700 pcs.     | 8715343034453 | 8715343034484 | 522160    |
|  | 3025565        | KLINION PROTECTION PE APRON 71 X 117 CM WHITE EMBOSSED 0.023 MM REF 522700-1       | 100 pcs.     | 500 pcs.     | 8715343066461 | 8715343066478 | 522700-1  |
|  | 3025566        | KLINION PROTECTION APRON 81 X 125 CM TRANSPARANT EMBOSSED 0.025 MM REF 522701-1    | 50 pcs.      | 500 pcs.     | 8715343066485 | 8715343066492 | 522701-1  |
|  | 3025739        | KLINION PROTECTION PE APRON 117 X 71 CM WHITE 18 MM EMBOSSED 0.018 MM REF 522700   | 100 pcs.     | 500 pcs.     | 8715343007846 | 8715343007990 | 522700    |
|  | 3025740        | KLINION PROTECTION APRON 125 X 81 CM TRANSPARANT EMBOSSED 0.0195 MM REF 522701     | 50 pcs.      | 500 pcs.     | 8715343007853 | 8715343008010 | 522701    |
|  | 3025741        | KLINION PROTECTION APRON 125 X 81 CM NATURAL EMBOSSED 0.0195 MM REF 522702         | 50 pcs.      | 500 pcs.     | 8715343007860 | 8715343008041 | 522702    |
|  | 3025742        | KLINION PROTECTION APRON 160 X 90 CM WHITE EMBOSSED 0.019 MM REF 522703            | 70 pcs.      | 420 pcs.     | 8715343007877 | 8715343008058 | 522703    |
|  | 3025743        | KLINION PROTECTION APRON 140 X 80 CM WHITE FOLDED 0.031 MM REF 522704              | 70 pcs.      | 420 pcs.     | 8715343007884 | 8715343008065 | 522704    |
|  | 3025744        | KLINION PROTECTION APRON ON ROLL 150 X 80 CM TRANSPARANT 0.0195 MM REF 522707      | 50 pcs.      | 500 pcs.     | 8715343445419 | 8715343445433 | 522707    |
|  | 3025745        | KLINION PROTECTION APRON FOLDED 150 X 80 CM TRANSPARANT 0.0195 MM REF 522708       | 40 pcs.      | 400 pcs.     | 8715343445426 | 8715343445440 | 522708    |
|  | 3049840        | KLINION PERSONAL PROTECTION DISPOSABLE LDPE APRON, WHITE 522709                    | 70 pcs.      | 420 pcs.     | 8715343065068 | 8715343065075 | 522709    |