

Product Identification

Product family	Warming blanket
Reference	4000002 (Former:85100)
Name	Blanket Full Body 195 x 100 cm
Sterilization	Non sterile
Brand	Sterisets®

Composition

Material	Material
Material	PP+PE nonwoven, Velcro
Tolerance	±2cm

Material

Test	Compliance	Result
PP + PE nonwoven	Internal supplier method	Yellow 50g/m ² , coated (yellow 30g PP + yellow 20g PE)
PP + PE nonwoven	Internal supplier method	White 50g/m ² , coated (white 30g PP + clear 20g PE)
Velcro tie	Internal supplier method	White, 1pc of 1x30cm
Length	Internal supplier method	195cm
Width	Internal supplier method	100cm
Airflow standard	Internal supplier method	5m/sec (+/-15%)
Weight	Internal supplier method	195g/pc (+/-5%)

Packaging

Packaging	Quantity	Type	UDI-DI (GTIN)
1 st packaging	1pc	Bag	08720598884685
2 nd packaging	10pcs	Inner box	18720598884682
3 rd packaging	40pcs	Carton box	28720598884689

Storage Conditions

Store in a dry place at room temperature. Do not expose to sunlight for long periods.

