

Product Identification	
Product family	Warming blanket
Reference	85200
Name	Blanket Lower Body 122 x 100 cm
Sterilization	Non sterile
Brand	Sterisets®

Composition	Material
Material	PP, PE and Velcro
Tolerance	±2cm

Material		
Test	Compliance	Result
PP + PE nonwoven	Internal supplier method	Yellow 50g/m ²
PP + PE nonwoven	Internal supplier method	White 50g/m ²
Velcro tie	Internal supplier method	White, 1cm wide x 30cm long
Length	Internal supplier method	122cm
Width	Internal supplier method	100cm
Airflow standard	Internal supplier method	4.81m/sec
Weight	Internal supplier method	125g/pc
Individual folding	Internal supplier method	22 x 26cm, yellow outside

Packaging			
Packaging	Quantity	Type	UDI-DI (GTIN)
1 st packaging	1pc	Bag	08720598881103
2 nd packaging	10pcs	Cardboard box	18720598881100
3 rd packaging	40pcs	Carton box	28720598881107

Storage Conditions

Store in a dry place at room temperature. Do not expose to sunlight for long periods.

