

Product Identification	
Product family	Warming blanket
Reference	4000004 (Former: 85310)
Name	Blanket Upper Body ARM OUT 192 x 57 cm
Sterilization	Non sterile
Brand	Sterisets®

Composition	Material
Material	PP+PE nonwoven, Velcro
Tolerance	±2cm

Material		
Test	Compliance	Result
PP + PE nonwoven	Internal supplier method	Yellow 50g/m ² , coated (yellow 30g PP + yellow 20g PE)
PP + PE nonwoven	Internal supplier method	White 50g/m ² , coated (white 30g PP + clear 20g PE)
Velcro tie	Internal supplier method	White, 4pcs of 1x30 + 2pcs of 2x2cm
Length	Internal supplier method	192cm
Width	Internal supplier method	57cm
Airflow standard	Internal supplier method	4.98m/sec (+/-15%)
Weight	Internal supplier method	125g/pc (+/-%)
Individual folding	Internal supplier method	22x26cm, yellow outside

Packaging			
Packaging	Quantity	Type	UDI-DI (GTIN)
1 st packaging	1pc	Bag	08720598884692
2 nd packaging	10pcs	Inner box	18720598884699
3 rd packaging	40pcs	Carton box	28720598884696

Storage Conditions

Store in a dry place at room temperature. Do not expose to sunlight for long periods.

