



## Safedraw®-P Safe Arterial Sampling Meritrans DTXPlus® Disposable Transducer Sets

Sterile, single use

PRODUCT CODE	NO. OF DISPOSABLE TRANSDUCERS	FLUSH/ FLOW RATE (mL/HR)	TOTAL PRESSURE LINE LENGTH (CM)	PRESSURE LINE(S) COLOUR	ADDITIONAL INFORMATION	SAFEDRAW SYRINGE VOLUME (mL)	SEPTUM MODEL	BOX/CASE CONTENTS (UNITS)
688413	1		183cm (153/30)	Clear	Single Admin Set. Single Transducer Set with Integrated Safedraw®-P Septum (Red)	2cc	1 x Inline	5/20
688410	1		213cm (183/30)	Red/Clear	Single Admin Set. Single Transducer Set with Integrated Safedraw®-P Septum (Red)	2cc	1 x Inline	5/25
688457	1		223cm (133/90)	Red/Clear	Single Admin Set. Single Transducer Set with Integrated Safedraw®-P Septum (Red)	3cc	1 x Inline	5/20
688538	1		78cm (63/15)	Red/Clear	Single Admin Set. Single Transducer Set with Integrated Safedraw®-P Septum (Red) with Velcro Strap for Arm mounting	1.5cc	1 x Inline	5/20
688403/ 689314*	1		A: 213cm (183/30) V: 150	A: Red/Clear V: Blue	Single Admin Set. Single Transducer Set with non-shunting stopcock and Integrated Safedraw®-P Septum (Red) (with Separately Packaged 150cm Blue pressure line)	2.5cc	1 x Inline	5/25
688407	1		A: 183cm (153/30) V: 125	A: Red/Clear V: Blue	Single Admin Set. Single Transducer Set with R.O.S.E., non-shunting stopcock, and Integrated Safedraw®-P Septum (Red) (with Separately Packaged 125cm Blue pressure line)	2cc	1 x Inline	5/20
688496	1		A: 193cm (163/30) V: 150	A: Red/Clear V: Blue	Single Admin Set. Single Transducer Set with non-shunting stopcock (with 150cm Blue pressure line) and Integrated Safedraw®-P Septum (Red)	2.5cc	1 x Inline (Red)	5/20
688416	2		A: 183cm (153/30) V: 150cm (125/25)	A: Clear V: Blue	Bifurcated Admin Set. Double Transducer Set with Integrated Safedraw®-P Septum (Red)	2cc	1 x Inline	5/20
688401	2		A: 213cm (183/30) V: 186cm (183/3)	A: Red/Clear V: Blue/Clear	Bifurcated Admin Set. Double Transducer Set with Integrated Safedraw®-P Septum (Red) and Integrated Safedraw®-P Stopcock Septum (Blue)	2.5cc	1 x Inline (Red) 1 x Stopcock (Blue)	5/25
688604	2		A: 223cm (133/90) V: 216cm (183/33)	A: Red/Clear V: Blue/Clear	Bifurcated Admin Set. Double Transducer Set with Integrated Safedraw®-P Septum (Red) and Integrated Safedraw®-P Stopcock Septum (Blue)	3cc	1 x Inline (Red) 1 x Stopcock (Blue)	5/25
688452	2		A: 223cm (133/90) V: 90cm	A: Red/Clear V: Blue	Trifurcated Admin Set. Double Transducer Set with Integrated Safedraw®-P Septum (Red) and Integrated Stopcock (Blue)	3cc	1 x Inline (Red)	5/20
688528	2		A: 223cm (133/90) V: 211cm (183/28)	A: Red/Clear V: Blue/Clear	Trifurcated Admin Set. Double Transducer Set with Integrated Safedraw®-P Septum (Red) and Integrated Safedraw®-P Stopcock Septum (Blue)	3.5cc	1 x Inline (Red) 1 x Stopcock (Blue)	5/20
688572	2		A: 233cm (143/90) V: 205cm (180/25)	A: Red/Clear V: Blue/Clear	Trifurcated Admin Set. Double Transducer Set with Integrated Safedraw®-P Septum (Red) and 2 x Sampling Stopcocks (Blue)	3cc	1 x Inline (Red)	5/20
688481	3		A: 213cm (183/30) V: 126cm (123/3) PA: 120cm	A: Red/Clear V: Blue/Clear PA: Yellow	Trifurcated Admin Set. Triple Transducer Set with Integrated Safedraw®-P Septum (Red) and Integrated Safedraw®-P Stopcock Septum (Blue)	3cc	1 x Inline (Red) 1 x Stopcock (Blue)	5/20
688396	1		193cm (163/30)	Red/Clear	Single Admin Set. Single 'Flexi' Transducer Set with non-shunting stopcock and Integrated Safedraw®-P Septum (Red)	2.5cc	1 x Inline	5/20
688389	1		193cm (163/30)	Red/Clear	Single Admin Set. Single 'Flexi' Transducer Set with Integrated Safedraw®-P Stopcock Septum (Red)	3cc	1 x Stopcock	5/20
688420/ 689311*	1		213cm (183/30)	Red/Clear	Single Admin Set. Single 'Flexi' Transducer Set with Integrated Safedraw®-P Septum (Red)	2.5cc	1 x Inline	5/25
688406	1		203cm (53/120/30)	Red/Clear	Single Admin Set. Single 'Flexi' Transducer Set with 2 x Integrated Safedraw®-P Septums (Red)	2.5cc	2 x Inline	5/20
688513	1		193cm (73/90/30)	Red/Clear	Single Admin Set. Single 'Flexi' Transducer Set with 2 x Integrated Safedraw®-P Stopcock Septums (Red)	3cc	2 x Stopcock	5/25
688408/ 689423*	1		143cm (128/15)	Clear	Single Paediatric Transducer Set with Integrated Safedraw®-P Septum (Red) <b>For use with Infusion Pump</b>	2.5cc	1 x Inline	5/25
688550/ 689312*	1		158cm (153/5)	Clear	Single Paediatric Transducer Set with Integrated Safedraw®-P Septum (Red) <b>For use with Infusion Pump</b>	2cc	1 x Inline	5/25
688600	1		158cm (153/5)	Clear	Single Paediatric Transducer Set with Integrated Safedraw®-P Septum (Red) <b>For use with Infusion Pump</b>	2cc	1 x Inline	5/25

\* Safedraw®-P DEHP-free sets