



Oxygen Supply Tubing

Oxygen Supply Tubing is a vital component of oxygen therapy, designed to ensure safe and efficient delivery of oxygen from the source to the patient.

Medline offers a variety of oxygen tubing in various lengths, softness, and connectors ensuring you have the right product to meet your clinical or financial goals.



Kink and crush-resistant tubing
This distinct tubing design features a five-channel inner lumen, which facilitates continuous oxygen flow.



Universal connectors
Convenient and economical because they attach directly to flowmeters, thus eliminating the need for nipple-and-nut adaptors.

Reference	Tubing Type		End connector		Pack factor
	Star Lumen	Length	Standard	Universal	
41113P*	X	210 cm	X		50
1118P*	X	420 cm	X		50
1119P*	X	760 cm	X		25
1120P*	X	1520 cm	X		25
1925P*	X	210 cm		X	50

*Di-2-ethylhexyl phthalate (DEHP)-free